

Uranyl Acetate Stain in Maleate Buffer (*en bloc*)

Warning: Uranyl acetate is toxic and radioactive. A discussion of the safe handling of this material can be found [here](#).

Materials required

Sodium Hydroxide (2N, Electron Microscopy Sciences (EMS), CAT #21170)

Maleate Acid (EMS, CAT # 18150)

Uranyl Acetate (EMS, CAT # 22400)

Distilled Water

Stock Solution

0.2N Sodium Hydrogen Maleate

For 250 ml:

1. Dissolve 5.8g of maleic acid in 50ml of 1N sodium hydroxide.
2. Dilute to 250ml with distilled water.

Procedure

For 800 ml:

1. To 107.6 ml 0.2 N sodium hydroxide
2. Add 200 ml 0.2 N NaH maleate (prepared above)
3. Dilute to near 800 ml with distilled water. Verify that the pH = 6.0 and adjust as necessary, then qc to 800ml.
4. Add 8g uranyl acetate to make a 1% stain solution.
5. Used to stain tissues after post-fixation by Osmium Tetroxide (described [here](#)).

Alternate Procedure

For 800ml:

1. Add 8g uranyl acetate to 800ml 0.2M Tris Maleate Buffer (pH 5.15, EMS, CAT # 11740).

Adapted from: Temple, J. (1929) Sodium Maleate-A Buffer for the P_H Region of 5.2 to 6.8. J. Am. Chem. Soc. **51**(6): 1754-1755.